

# **RNA injection protocol for *Xenopus tropicalis***

**09/18/2018**

**Hajime Ogino, Takeshi Igawa, and Nanoka Suzuki**  
**(Amphibian Research Center, Hiroshima University)**

**Haruki Ochi**  
**(Faculty of Medicine, Yamagata University)**

---

---

## EQUIPMENT

---

### Pressure microinjector

- HARVARD APPARATUS PLI-100 Pico-Injector
- NARISHIGE IM-300 Microinjector
- NARISHIGE IM-31 Microinjector

IM-31 is the cheapest injector, but appears to be sufficient for RNA injection.

---

### Puller

- NARISHIGE PN-31
- 

### Manipulator

- NARISHIGE M-152
- NARISHIGE MN-151
- NARISHIGE YOU-1


Choose a favorite one.

---

### Injection dish

- Non-toxic Modeling Clay (Van Aken Plastalina)
- 6 cm dish (VALMARK 2901)

Put a layer of the modeling clay in the bottom of a petri dish. Make grooves or holes on its surface to hold fertilized eggs.




### Culture dish

- 6-well dish or a petri dish (100 mm × 15 mm)
- 0.7 – 1 % Agarose in 0.1 x MMR (pH7.5)

Coat the dish bottom with agarose to prevent sticking of embryos.

This agarose-coating is not necessary if 0.1% BSA is included in the culture medium.



---

### **Injection needles**

- Borosilicate Glass Capillaries (World Precision Instruments (WPI), TW100-4)
- NARISHIGE PN-31 Puller

#### Puller setting condition

Heater	80
Magnet sub	30
Magnet main	90-100

---

### **Injection needle calibration**

- Cut the needle tip using a sterile forceps.
- Fill 1  $\mu$ l of RNase-free water into the needle.
- Set the injection time and pressure (1 second and 10 psi).
- Eject all of the RNase-free water. If it takes about 40 – 60 seconds, you can use that needle for injection into *X. tropicalis* embryos.

---

### **Plastic pestle for 1.5 ml microcentrifuge tube** (Kontes, 749521-1590)

**Air incubator** (temperature should be adjustable between 12 –25°C)

---

---

## **ANIMALS**

---

Researchers can purchase *X. tropicalis* from the “National Bioresource Project of Japan” ([http://home.hiroshima-u.ac.jp/amphibia/xenobiores\\_en/iweb\\_en/Top.html](http://home.hiroshima-u.ac.jp/amphibia/xenobiores_en/iweb_en/Top.html)). The four inbred wild-type strains, Nigerian A, Golden, Nigerian H, and Ivory Coast, are available. The Golden strain is very tough and appear to be most suitable for routine injection experiments

---

---

## REAGENTS

---

<b>10 x MMR</b> (1 L)	<u>from stock solutions</u>	<u>in grams</u>
1 M NaCl	250 ml 4 M NaCl	58.44 g NaCl
20 mM KCl	20 ml 1 M KCl	1.49 g KCl
10 mM MgCl <sub>2</sub>	10 ml 1 M MgCl <sub>2</sub>	2.03 g MgCl <sub>2</sub> •6H <sub>2</sub> O
50 mM HEPES	50 ml 1 M HEPES	11.92 g HEPES
20 mM CaCl <sub>2</sub>	20 ml 1 M CaCl <sub>2</sub>	2.94 g CaCl <sub>2</sub> •2H <sub>2</sub> O

Add distilled water to 1 Liter. Adjust pH to 7.5 by adding approximately 4 ml of 5M NaOH (Nacalai 95539-85). Sterilize by autoclave or filtration if made from dry stocks. Otherwise just combine sterile solutions + sterile H<sub>2</sub>O in a sterile bottle. Store at RT.

---

### **0.1 x MMR (pH7.5)**

- 990 ml Distilled water
- 10 ml 10 x MMR

Re-adjust the pH to 7.5 by adding about 10 ul of 5M NaOH.

---

### **0.1 x MMR (pH7.5) + Gentamicin (embryo culture medium)**

- 990 ml Distilled water
- 10 ml 10 x MMR
- 1 ml 50 mg/ml Gentamicin Sulfate (Biowest P4020-5GR)
- 100 µl 0.5 % Phenol Red solution (Sigma P0290)

Re-adjust the pH to 7.5 by adding about 20 ul of 5M NaOH (gentamicin decreases the pH).

---

### **Dejelling solution (2% L - cysteine in 0.1 x MMR, pH 8.0)**

Prepare freshly on the day of injection.

- 50 ml 0.1 x MMR
- 1 g L – cysteine (Wako 326-0615)

Adjust the pH to 8.0 by adding 600 – 700 µl of 5 M NaOH. Such higher pH is crucial for dejelling.

---

**Injection medium (6% Ficoll 400 and 0.1% BSA in 0.1 x MMR, pH 7.5)**

- 500 ml 0.1 x MMR, pH7.5
- 30 g Ficoll 400 (GE Healthcare 17-0300-50)
- 5 ml 10% BSA (pH7.5)

Dissolve Ficoll 400 powder slowly (gently swirl one overnight).

The pH may need to be re-adjusted to 7.5.

Make 50 ml aliquots, and store at -20°C.

---

**Anesthetic solution (0.2% MS222 in 0.1 x MMR)**

- 990 ml Distilled water
- 10 ml 10 x MMR
- 2 g MS222 (Ethyl 3-aminobenzoate methanesulfonate, Sigma E10521)

Re-adjust the pH if necessary.

---

**Testis storage buffer**

Leibovitz L15 medium is convenient. One may keep *X. tropicalis* testes for up to two days at 14°C.

Preparation of the 190% Steinberg solution takes time, but it may be efficient for longer storage.

**1. Leibovitz L15 medium (Sigma L4336)**

Follow the manufactures instructions for preparation.

Supplementing with gentamycin (Final 50 ug / ml) may be effective for longer storage of the testes.

**2. 190% Steinberg solution (pH 7.8)**

- 110 mM NaCl
- 1.29 mM KCl
- 0.64 mM  $\text{Ca}(\text{NO}_3)_2 \cdot 4\text{H}_2\text{O}$
- 1.58 mM  $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$
- 8.78 mM Tris-base

Add gentamycin (Final 50 ug / ml) and 1/10,000 vol. of 0.5 % Phenol Red solution.

Adjust the pH to 7.8, then sterilize by filtration. Store at 4°C.

---

**10% BSA (pH7.5)**

- 50g BSA (Fraction V, Roche 735094)
- 400 ml Distilled water

Adjust the pH to 7.5, fill up to 500 ml.

Sterilize by filtration, make 50 ml aliquots, and store at -20°C.

---

**Testis homogenizing buffer (0.1% BSA in 1 x MMR)**

- 45 ml Distilled water
- 5 ml 10 x MMR
- 0.5 ml 10% BSA

Make 1 ml aliquots, and store at -20°C.

---

**Capped RNAs for injection**

One may prepare Capped RNAs by in vitro transcription using a mMMESSAGE mMACHINE kit (Thermo Fisher Scientific) or equivalent kits/methods. For details of the RNA preparation method, please see published protocols (Early Development of *Xenopus laevis*, A Laboratory Manual, CSHL press, ISBN 978-087969942-0). The plasmid template may be constructed by introducing your experimental cDNA into the MCS of pCS2+ vector, in which the MCS is flanked by a SP6 promoter and a SV40 poly-A signal sequence (Rupp, R.A., et al., Genes Dev. 8, 1311–1323, 1994).

---

**HCG (Human chorionic gonadotropin)** (Aska Pharmaceutical 45254000, or equivalent)

---

---

## PROCEDURE

- 1 In the afternoon of the day before the RNA injection, prime *X. tropicalis* females about 20 hr prior to egg laying with 15 units of HCG (15 units/100 ul/frog), and keep them at 25°C in holding tanks.
- 2 On the morning of the day of RNA injection, boost the females with 100 units of HCG approximately 4 hr prior to the desired time of egg laying. Equilibrate necessary reagents and culture dishes at room temperature (22-25°C).
- 3 Euthanize *X. tropicalis* males to isolate testes. Pre-warm the anesthetic solution to room temperature, and submerge the males in this solution for between 30 minutes and 1 hour. After euthanasia, dissect out the testes and transfer them to Leibobitz L-15 medium or 190% Steinberg solution. For details of the testis dissection, one should follow standard procedures as described (Early Development of *Xenopus laevis*, A Laboratory Manual, CSHL press, ISBN 978-087969942-0).

## CAUTION

Please follow animal protocols of your institute for frog handling.

- 4 When the female frogs start laying eggs in their containers, start *in vitro* fertilization. Homogenize a testis with 0.5-0.75 ml of the testis homogenizing buffer in a 1.5 ml microcentrifuge tube using a plastic pestle. Gently squeeze 500-600 eggs out of the female in a dry Petri dish, remove all water from the eggs, and then mix them with a part of homogenized testis using a wide-pore pipet. After two minutes, flood eggs with 50 ml of 0.1×MMR (room temperature) in the Petri dish and further incubate them at room temperature.
- 6 If the fertilization is successful, contraction of the pigmented animal cap occurs after 10-15 minutes. Thirty minutes after the addition of homogenized testis, replace the flooded 0.1×MMR with 50-100 ml of the dejelling solution and swirl eggs until their jelly coat is completely removed (2-5 minutes). After dejelling, wash the dejellied eggs several times with 50-100 ml of 0.1×MMR.
- 7 Fill injection dishes with the injection medium, and lay the required number of fertilized eggs in the grooves or holes on the dish bottom.
- 8 Load 1-2 ul of the capped RNA solution into an injection needle. At 1-cell stage, one can inject up to 4 nl per embryo. At 2-cell stage, one can inject up to 2 nl per one blastomere of an embryo. For details

of the microinjection itself, one should follow standard *Xenopus* procedures as described (Early Development of *Xenopus laevis*, A Laboratory Manual, CSHL press, ISBN 978-087969942-0).

- 9** After the injection, incubate the embryos at 22-25 °C in the injection medium for about 2 hours to let healing.
  
- 10** Transfer embryos from the injection dish to Petri dishes filled with 0.1 x MMR (pH7.5) containing gentamycin, and further incubate them at 22-25 °C. Do not place too many embryos in the same dish (up to 100 embryos per Petri dish (100 mm × 15 mm)).
  
- 11** Next morning, the embryos will reach stages 20-25. To keep their development, change the culture medium every day, and remove dead embryos.